VDRL Cardioliipin Antigen

Intended Use

VDRL (Venereal Disease Research Laboratory) Cardioliipin Antigen is a precipitation reagent that aids in the detection and quantification of lipoidal antibodies in serum.

The test was developed for testing individual samples, not for pooled samples.

The product is for *in vitro* diagnostic use only.

Summary and Explanation

The detection and quantification of lipoidal antibodies is of diagnostic importance in connection with the diagnosis of syphilis, primarily for follow-up and assessment of the effectiveness of therapy¹⁻⁵. The findings obtained must be considered in conjunction with the results of treponema-specific tests (Enzygnost® Syphilis, IgM-FTA-ABS, FTA-ABS, indirect hemagglutination).

Refer to the relevant literature for the evaluation of findings.

Principle of the Method

Serum or liquor is mixed with cardioliipin antigen on a slide. The presence of anti-lipoidal antibodies causes agglutination (formation of precipitate) with cardioliipin antigen.

If no anti-lipoidal antibodies are present, the suspension remains cloudy throughout.

Reagents

Materials provided in the VDRL kit:

VDRL Cardioliipin Antigen	1 pcs.
Dilution liquid	1 pcs.
Instructions for use	1 pcs.

Composition

VDRL Cardioliipin Antigen, bovine, in alcoholic solution:

Cardioliipin	0.3 g/L
Lecithin	2.1 g/L
Cholesterol	9.0 g/L

Dilution fluid

Formaldehyde	0.58 g/L
Na ₂ HPO ₄ x 12 H ₂ O	93 mg/L
KH ₂ PO ₄	170 mg/L
NaCl	10.0 g/L

Warnings and Precautions

- For *in vitro* diagnostic use only.
- VDRL Cardioliipin is contained in an alcoholic solution and is highly flammable. The container must be kept tightly closed and removed from sources of ignition.
- It is advisable to wear protective gloves throughout the entire test procedure.
- For the disposal of solid infectious materials, autoclaving for at least 1 hour at +121 °C is recommended. Disinfectants suitable for inactivating human pathogens may also be used. The concentrations and times specified by the manufacturer must be observed.

Preparations of Reagents

Preparation of the Antigen Suspension

Put 0.4 mL dilution fluid into a glass vial and then slowly add 0.5 mL of antigen while agitating.

Do not allow the pipette to come into contact with the dilution fluid.

Squeeze out the last few drops and then gently agitate the fluid by swaying it from side to side a few seconds more.

Add 4.1 mL dilution fluid, close the vial, and shake vigorously.

The suspension must be used within one day.

Preparation of serum

Inactivate sera for 30 minutes at +56 °C ±2 °C. Samples that have already been inactivated must be inactivated again after more than 4 hours for 10 minutes at +56 °C ±2 °C.

Liquors are not inactivated.

Storage and Stability

All components of the test kit may be used up to the date indicated on the label if stored at a temperature between +15 to +25 °C.

For complete stability and storage data for opened reagents or reagents diluted for use, see the following table.

Reagent	State	Storage	Stability
VDRL Cardioliipin Antigen	after opening tightly closed	+15 to +25 °C	up to expiry date
VDRL-Cardioliipin-Antigen	diluted ready for use	+15 to +25 °C	1 day
Dilution liquid	after opening tightly closed	+15 to +25 °C	up to expiry date

Materials required but not provided

Lues Standard Serum (Code No. ORDN) as a control for the reaction and for quantitative evaluation of cardioliipin antibodies.

Slides with wells.

Specimens

Suitable specimens for testing are individual serum samples obtained by standard laboratory techniques.

The samples should be stored for no longer than 3 days at +2 to +8 °C. If the samples are to be stored for a longer period of time, they must be frozen.

Inactivate the sera before testing for 30 minutes at +56 °C ±2 °C.

Procedure

Slide Test with Serum

- Dispense and distribute 50 µL inactivated serum onto a slide with wells (14 mm).
- Add 17 µL antigen suspension with a piston-type pipette.
- Agitate for 4 minutes on a rotation device (180 ± 50 rpm and 1.9 cm circular movement).

The assessment is performed immediately upon finishing the rotation with a microscope at 100 fold magnification:

no agglutination	non-reacting
agglutination not clearly recognizable	questionable reaction
small but clear agglutination	weak reaction
strong agglutination	reacting

Internal Quality Control

Lues Standard Serum is used as a control for internal quality control.

Quantitative determination of antibody content

For the determination of cardioliipin antibodies, the slide test is performed with 50 µL each of a geometric dilution series of the inactivated serum. Dilution fluid is recommended for the dilution of serum.

The corrected mean titer or the number of International Units per mL is based on comparative titration of the Lues Standard Serum according to the recommendations of the German Federal Health Authority. The titer is the reciprocal value of the last dilution level to be determined as "reacting".

Limitations of the Procedure

- Incompletely coagulated sera and microbially contaminated test samples should not be used. Any particulates (e.g. fibrin clots, erythrocytes) contained in the samples should be removed prior to assay.
- If thawed samples are used, ensure that the material is thoroughly homogenized.
- The reagent is only suitable for testing individual samples, not for pooled samples.
- The instructions given for the preparation of the antigen suspension must be followed very closely.

Specific Performance Characteristics

VDRL Cardioliipin Antigen is suitable for the detection of anti-lipoidal IgG and IgM antibodies.

The performance of each lot of reagent is tested with clinically known positive and negative controls.

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Enzygnost is a registered trademark of Dade Behring Marburg GmbH in the USA, Germany and other countries.

 Dade Behring Marburg GmbH
Emil-von-Behring-Str. 76
D-35041 Marburg
www.dadebehring.com

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